







IPL SYSTEMS QUOTATION

MODEL	DESCRIPTION	PACKING INFORMATION	PRICE(AUD)	IMAGE
EROSE-YA	<p>Spectrum range: 420-1200nm; 520-1200nm; 620-1200nm Option: 560-1200nm; 690-1200nm; 755-1200nm Energy density: 10-50J/cm² Pulse sequence: 2,3,5 pulse Pulse intermission: 5-60ms Pulse width: 2-15ms Repeat frequency: 0.3-1Hz Delivery: Sapphire Coupling Light spot: 8×40mm² Cooling system: Water cooling, forced-air cooling and semiconductor cooling</p> <p>Cooling temperature: -4 - 0℃ Power supply: 220V / 50Hz or 110v / 60Hz</p>	<p>Wooden box Packing N.W:58kgs; G.W:96kgs 74*60*134cm/Wooden box; 63*46*115/Machine size</p>	14,500	
EROSE-YC	<p>Spectrum range: 420-1200nm; 520-1200nm; 620-1200nm Option: 560-1200nm; 690-1200nm; 755-1200nm Energy density: 10-50J/cm² Pulse sequence: 2,3,5 pulse Pulse intermission: 5-60ms Pulse width: 2-15ms Repeat frequency: 0.3-1Hz Delivery: Sapphire Coupling Light spot: 8×40mm² Cooling system: Water cooling, forced-air cooling and semiconductor cooling</p> <p>Cooling temperature: -4 - 0℃ Power supply: 220V / 50Hz or 110v / 60Hz</p>	<p>Wooden box Packing N.W:57kgs; G.W:93kgs 71*74*130cm/Wooden box; 63*68*118/Machine size</p>	13,000	

IPL SYSTEMS QUOTATION

MODEL	DESCRIPTION	PACKING INFORMATION	PRICE(AUD)	IMAGE
EROSE-Y3	<p>Spectrum range: 520-1200nm; 620-1200nm Option: 420-1200nm; 560-1200nm; 690-1200nm; 755-1200nm Energy density: 10-50J/cm² Pulse sequence: 1,2,3,4,5 pulse Pulse intermission: 5-60ms Pulse width: 2-15ms Repeat frequency: 0.3-1Hz Delivery: Sapphire Coupling Light spot: 8×40mm² Cooling system: Water cooling, forced-air cooling and semiconductor cooling Cooling temperature: -4 - 0℃ Power supply: 220V / 50Hz or 110v / 60Hz</p>	<p>Wooden box Packing N.W:44kgs; G.W:64kgs 61*58*123cm/Wooden box; 50*50*98/Machine size 1.accessorie ctn G.W:7kgs 55x25x33cm/ctn</p>	10,200	
HM-IPL-B3	<p>Spectrum range: 420-1200nm; 510-1200nm; 560-1200nm;640-1200nm Option: 620-1200nm;755nm-1200nm Option: 585-1200nm; 690-1200nm Pulse sequence: 1,2,3,4,5 pulse Pulse intermission: 5-60ms Pulse width: 2-15ms Repeat frequency: 0.3-1Hz Delivery: Sapphire Coupling Power: 1200W Light spot: 8×40mm² Cooling system: Water cooling, forced-air cooling and semiconductor cooling Cooling temperature: -4 - 0℃ Power supply: 220V / 50Hz or 110V 60Hz</p>	<p>Aluminium Packing N.W:31kg G.W:42kg 68*59*51/Aluminum box 58*47*34/Machine size</p>	8,000	

IPL SYSTEMS QUOTATION

MODEL	DESCRIPTION	PACKING INFORMATION	PRICE(AUD)	IMAGE
HM-IPL-B2	<p>Spectrum range: 420-1200nm; 510-1200nm; 560-1200nm;640-1200nm Option: 620-1200nm;755nm-1200nm Option: 585-1200nm; 690-1200nm Pulse sequence: 1,2,3,4,5 pulse Pulse intermission: 5-60ms Pulse width: 2-15ms Repeat frequency: 0.3-1Hz Delivery: Sapphire Coupling Power: 1000W Light spot: 8×40mm² Cooling system: Water cooling, forced-air cooling and semiconductor cooling Cooling temperature: -4 - 0℃ Power supply: 220V / 50Hz or 110V 60Hz Spectrum range: 420-1200nm; 510-1200nm; 560-1200nm;640-1200nm Option: 585-1200nm; 690-1200nm Energy density: 10-50J/cm² Pulse sequence: 1,2,3,4,5 pulse Pulse intermission: 5-60ms Pulse width: 2-15ms Repeat frequency: 0.3-1Hz Delivery: Sapphire Coupling Power: 800W Light spot: 8×40mm² Cooling system: Water cooling, forced-air cooling and semiconductor cooling Cooling temperature: -4 - 0℃ Power supply: 220V / 50Hz or 110V 60Hz</p>	<p>Aluminium Packing N.W:31kg G.W:42kg 68*59*51/Aluminum box 58*47*34/Machine size</p>	7,500	
HM-IPL-B1	<p>Spectrum range: 420-1200nm; 510-1200nm; 560-1200nm;640-1200nm Option: 585-1200nm; 690-1200nm Energy density: 10-50J/cm² Pulse sequence: 1,2,3,4,5 pulse Pulse intermission: 5-60ms Pulse width: 2-15ms Repeat frequency: 0.3-1Hz Delivery: Sapphire Coupling Power: 800W Light spot: 8×40mm² Cooling system: Water cooling, forced-air cooling and semiconductor cooling Cooling temperature: -4 - 0℃ Power supply: 220V / 50Hz or 110V 60Hz</p>	<p>Aluminium Packing N.W:21kg G.W:33kg 59*45*47/Aluminum box 47*37*34/Machine size</p>	7,000	

IPL SYSTEMS QUOTATION

Accessory: each of system contains: Manual, Video, 1 Protective glasses for Doctor, 1 Eye patch for Patient, 1 Pipe, 1 Funnel, and 6 Bottles of GEL.

REMARK

1	Currency	Price is quoted in AUD are not inclusive of GST
2	Price Term	The prices quoted above are current as at 1 st January 2008 and may be subject to change
3	Delivery Time	Within 10 – 15 days after receipt of payment
4	Payment Term	Direct Bank Transfer in advance
5	Certificate	SFDA; CE
6	Warranty	Standard warranty of system unit is 12 months ex factory

IPL APPLICATION

- 1) Hair removal
- 2) Skin rejuvenation
- 3) Pigmentation treatment
- 4) Eliminate wrinkles
- 5) Acne clearance
- 6) Vascular treatment

IPL Systems

24 Fremantle Rd
Port Noarlunga South,
South Australia 5167
Phone: 618 8327 0845 • Fax: 618 8327 0845

Business Hours:

Monday - Friday: 9:00 am - 5:00 pm CST

Contact Email Address:

info@ipl-systems.com